

4-23-2011

FT-950 MIC SETTING FOR HEIL GM-4 GM-5

91	200	-100
92	-8	-5
93	2	1
94	900	800
95	-2	5
96	10	2
97	2600	3200
98	8	8
99	3	1

PROSSIOR SETTING

100	200	200	100
101	-5	-8	-5
102	4	2	1
103	1500	900	800
104	3	-6	5
105	5	2	2
106	2800	2100	3200
107	5	-8	8
108	3	2	1
109	5	5	7

The best way for you to come up with what sounds good is to use a headset. Under the power switch on/off push the button labeled MONI .

Listen to how you sound, next push the PROC button in once listen again if it doesn't sound the greatest adjust 91 through 99 from the menu for the best audio.

The processor is pretty touchy